

CIRCAID Custom Armsleeve & Glove Form

Fax Orders: 888-840-0939 Email: customs@mediusa.com

☐ Exact Reorder Order

Number: _____

RMA #: _____

medi

Customer Name _____

Customer No. _____ Purchase Order No. _____

Patient Name _____ Fax _____

Billing Address _____

Shipping Address _____

Telephone _____ Measured by _____

Order Date _____ Email _____

Shipping method: ☐ Standard (max. 5 days if complete order is received by 3:00pm EST) ☐ Second Day (extra charge) ☐ Next Day (extra charge)

Credit Card Info _____

juxtafit® essentials arm				Color ☐ Beige ☐ Black NEW			
	Qty Left	Qty Right	Standard arm sleeve (C-G)	Partial sleeve (C-E)	Partial sleeve (E-G)	Partial sleeve (C-F)	Made of extra-firm material
armsleeve							

juxtafit® essentials glove				Color ☐ Beige ☐ Black NEW			
*open palm is recommended for palm circumferences greater than 24cm	Qty Left	Qty Right	Fingerless glove	Extend wrist band			
glove				☐ 3cm	☐ 6cm		
open palm glove				☐ 3cm	☐ 6cm		
glove w/ dorsum strap				☐ 3cm	☐ 6cm		
open palm glove w/ dorsum strap				☐ 3cm	☐ 6cm		
hand wrap							

Cover up color

(indicate qty.)

____ black

(default)

____ beige

____ pink

Additional Options

Undersleeves (indicate qty.)

____ standard

____ silver

____ wide

*New custom options and add-ons may incur additional costs.

graduate™						
	Qty Left	Qty Right	Foam lateral rise	Stop foam at C	Tack thumb hole	Band locks
armsleeve			☐ 2.5cm ☐ 5cm			
graduate™ foam liner color options						
Interior			Exterior			
☐ Beige (default)			☐ Black (default) ☐ Beige			

Norton Hand Solution	Qty Left	Qty Right
shell and black glove		
shell and black glove with open palm		

Notes on taking measurements for custom-made circaid inelastic products

- Measurements for compression garments should not be taken until the best possible decongestion effort has been achieved. Circumference and length measurements are taken in a distal to proximal sequence. Measurements should be taken with the arm extended straight, such as with a patient lying flat on a table or standing.
- It is essential to mark the measuring points on the arm so that the circumference and length measurements are taken at the same point.
- The circumference measurements are taken without any tension as the products are adjustable and will accommodate some changes in size. Measurements are to be taken with skin measurements. Skin measurements should be taken loosely without tension.

Please take measurements without tension!


Patient Name:

Date: _____

Arm measurement

Left in cm.		Right in cm.
	55	
	50	
	45	
	40	
	35	
	30	
	25	
	20	
	15	
	10	
	05	
	C1	
	C wrist	
	B Palm	

Please hold the
arm straight.




Measurements must be every 5cm from the starting point (C).

Arm measurement Lengths

Left in cm.		Right in cm.
	ℓ C-G	
	ℓ E-G	
	ℓ C-F	
	ℓ C-E	



*Only C-E and C-G required for standard arm sleeve. Other measurements only needed for partial arm sleeves depending on sleeve end point.

Please measure the length on the inside of the arm in a straight line.





Hand measurement Circumferences

Left in cm.		Right in cm.
	1ez	
	2ez	
	3ez	
	4ez	
	5ez	
	1ex	
	2ex	
	3ex	
	4ex	
	5ex	
	cA	
	cB Palm	
	cC Wrist	
	cC1*	
	10*	



Hand measurement Lengths

Left in cm.		Right in cm.
	ℓ_1	
	ℓ_2	
	ℓ_3	
	ℓ_4	
	ℓ_5	
	ℓ_6	
	ℓ_7	
	ℓ_8	
	ℓ_9^1	
	ℓ_{10}^1	



*C1 and 10 required for extended wrist band on glove. 10 = 10 cm above C (wrist) measurement.

¹When measuring these points, press the palm flat on top of the measuring tape on a hard flat surface and measure the widths from the lateral edges of the hand.